

REDEExtract-N-Amp(TM) PCR ReadyMix(TM), Ready-to-use 2X PCR Master Mix with Loading Dye

REDEExtract-N-Amp(TM) PCR ReadyMix(TM),
Ready-to-use 2X PCR Master Mix with Loading
Dye

Identificativi Prodotto

Reference: R4775-1.2ML



Descrizione

Dati Tecnici

Purity Grade: Ready-to-use 2X PCR Master Mix with Loading Dye
Shipped on ICE/no ICE: WICE
SAP Pack with Pack Qty: 1X1.2ML
Storage Temperature: -20C