


SYBR(R) Green JumpStart(TM) Taq ReadyMix(TM), for quantitative PCR, MgCl₂ in buffer

SYBR(R) Green JumpStart(TM) Taq
ReadyMix(TM), for quantitative PCR, MgCl₂ in
buffer

Identificativi Prodotto

Reference: S4438-100RXN



Descrizione

Dati Tecnici

Purity Grade: for quantitative PCR, MgCl₂ in buffer

Shipped on ICE/no ICE: WICE

SAP Pack with Pack Qty: 1X100RXN

Storage Temperature: -20C