

Extract-N-Amp(TM) Tissue PCR Kit, sufficient for 1000 extractions, sufficient for 1000 amplifications

Extract-N-Amp(TM) Tissue PCR Kit, sufficient for 1000 extractions, sufficient for 1000 amplifications

Identificativi Prodotto

Reference: XNAT2R-1KT



Descrizione

Dati Tecnici

Purity Grade: sufficient for 1000 extractions, sufficient for 1000 amplifications

Shipped on ICE/no ICE: WICE

SAP Pack with Pack Qty: 1X1KIT

Storage Temperature: -20C